


New product release

Quartz wool discs suitable for elemental analysers

- ü Quick and simple reaction tube packing
- ü No waste
- ü Neat, symmetrical separation of reagents
- ü Available in three convenient sizes
- ü Reduced handling problems
- ü High purity

Discs are available in three different diameters – 15mm, 16.5mm and 22mm to fit most popular ‘micro’ elemental analysers currently in use. Thickness is 5mm approx. No more messing about with

trying to force small tufts of quartz wool into reaction tubes – simply insert the disc and push down onto the reagent to form a quick and effective separation. Handling problems from dust and fibres are reduced dramatically from using the traditional quartz wool.



Ordering details

Cat No	Description	Pack
B1250	Quartz Wool Discs 15mm	25
B1251	Quartz Wool Discs 16.5mm	25
B1252	Quartz Wool Discs 22mm	25

Recommended sizes

Manufacturer	Model	Cat No	Diameter
Costech	All	B1250	15mm
Europa	All	B1251	16.5mm
Eurovector	All	B1250	15mm
Perkin Elmer	2400	B1252	22mm
	2410	B1252	22mm
Thermo/Fisons/Carlo Erba	1108	B1250	15mm
	1110	B1250	15mm
	1112 (CHNOS)	B1250	15mm
	1112 (N/Protein)	B1252	22mm
	1500 series 1	B1251	16.5mm
	1500 series 2	B1250	15mm
	2000	B1250	15mm
	2100	B1250	15mm
	2500	B1250	15mm



Elemental Microanalysis
1 Hameldown Road
Okehampton
Devon, EX20 1UB, UK
+44 (0)1837 54446
info@microanalysis.co.uk
www.microanalysis.co.uk

Nov 2005

Or your local dealer: