

# New Product Release



From Left to Right:  
D2301, D2309, D2311.

## Lightweight Silver Capsules

- Reduce ash accumulation
- Great economy alternative
- Variety of Heights and Diameters

## AVAILABLE FROM:

Elemental Microanalysis Ltd

1 Hameldown Road Okehampton Devon  
EX20 1UB UK

T +44 (0)1837 54446

F: +44 (0)1837 54544

info@microanalysis.co.uk

[www.microanalysis.co.uk](http://www.microanalysis.co.uk)

[www.elementallab.co.uk](http://www.elementallab.co.uk)

Or your local dealer:

Elemental Microanalysis Ltd is pleased to introduce **NEW Lightweight Silver Capsules**

These new silver capsules have two significant advantages:

- Reduced accumulation of ash in the reaction tube, allowing more samples to be analysed before ash removal
- Significantly lower prices than the 'standard weight' versions – doubly important when silver has increased so much in cost in the last few years.

Products available  
and prices:

Part No	Size (Height X Diameter)	Pack	Weight (mg)	Volume ( $\mu$ l)
D2300	4 x 3.2mm	500	7.5	32
D2301	4 x 3.2mm	250	7.5	32
D2302	5 x 3.5mm	500	11.7	48
D2303	5 x 3.5mm	100	11.7	48
D2306	6 x 4mm	500	16.5	75
D2307	6 x 4mm	100	16.5	75
D2308	8 x 5mm	250	36	157
D2309	8 x 5mm	100	36	157
D2310	12.5 x 5mm	100	72.8	245
D2311	11 x 8mm	100	73.3	553
D2312	10.5 x 9mm	100	73.3	668