



high quality consumables for Trace elemental analysis



Product Sheet ~ Shimadzu TOC

Chemical Reagents

Cat #	Description	OEM #
B1605	High Sensitivity Catalyst 630-00996	630-00996 630-00996-00
B1606	High Sensitivity TC catalyst set 638-92070-01	638-92070-01
B1600	IC reagent 250ml 630-00710	630-00710
B1246	Platinum Catalyst 0.5% Pt 30gm	017-42801-00 017-42801-01 638-60116 638-60116-00
B1196	Platinum Catalyst 5% Pt 03002262 30gm	-
B1601	Potassium hydrogen phthalate 630-635-01 25g	630-635-01
B1017	Silica (Quartz) Wool Very Fine 5gm	630-00557 630-00557-00
B1602	Sodium carbonate 630-00962-01 25g	630-00962-01
B1603	Sodium hydrogen carbonate 630-00963-01 25g	630-00963-01
B1607	TC Catalyst Set 638-92069-02	638-92069-02

O-rings and operational accessories

Cat #	Description	OEM #
E7074	4-Way Valve 638-56079	638-56079-00
E7070	5ml Sample vials Round bottom no cap Shimadzu pack of 100	-
E7071	5ml Sample vials Round bottom with cap and seal Shimadzu pack of 100	-
E7054	5ml syringe for Shimadzu TOC 638-59161-10	638-59161-10
E7072	Chart Paper 108x50m 630-089-13-01 Pack of 10	200-91532-02 630-08913-01
E1550	O-Ring 4D P10A Shimadzu 036-11209	036-11209-84
E1551	O-Ring 4D P22 Shimadzu 036-11222	036-11222-84
E1552	O-Ring AS 568A-116 Shimadzu 630-01566	630-01566-00
E1553	O-Ring Teflon P10 TOC5000 036-11408	036-11408-84
D7018	Platinum Gauze Disc 630-00105-01 pack of 2	630-00105-01
E2158	Silicone Grease 501-241 100gm	-
E1554	Sleeve set 6F-T Front and back ferrule 035-62994-03	035-62994-03
E7073	Slider for TOC 631-32003	631-32003-00
E7051	Syringe 2500ul 630-02591-07	228-25237-06 630-02591-07
E7050	Syringe 250ul 630-02591-04	220-95277-00 630-02591-04
E7053	Teflon Syringe Tips for 2500ul Syringe Pack of 5 638-60734-08	638-60734-08
E7052	Teflon Syringe Tips for 250ul Syringe Pack of 5 638-60734-05	638-60734-05

Prepacked reactors and consumables kits

Cat #	Description	OEM #
		630-00999
X4010	CO2 Scrubber Shimadzu TOC 630-00999	630-00999-00
X4000	Halogen Scrubber Shimadzu TOC 630-00992	630-00992 630-00992-00

Quartz, Glassware and Furnace tubes

Cat #	Description	OEM #
C8060	Combustion Tube Shimadzu TOC 631-41524	631-41524-00 638-41323 638-41323-00
C8061	Shimadzu combustion tube 638-42076	638-42076

Placing your order: We are happy to accept your order by email, telephone, fax or post. We welcome all orders, large or small, with no minimum order requirement. Visa and Mastercard payments are accepted or we can open an account for you – please enquire. Most orders are despatched from stock within 24 hours.

Quality: All products are manufactured in our ISO9001:2000 approved facility and are covered by our 'Total Satisfaction Guarantee'.

Can't see what you need? We offer thousands of products for all elemental analysers. If you can't see what you need here the chances are we will have the product available. Just contact our friendly sales team who will track down what you need.

OEM Part numbers: The use of OEM part numbers is intended for convenience only and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and suitable for use in the stated instrument. All trade marks acknowledged.