

## New Product Release

### Combustion Tube Conversion Kit for Leco® CS844/CS744

Elemental Microanalysis offers a conversion kit, [E4356](#), that adapts the Leco® CS844/CS744 analysers to be able to use a much lower cost combustion tube (C4356) than the Leco® equivalent tube.

Customers need only purchase this easy-to-use kit once, to adapt the analyser, and then the combustion tube [C4356](#) from Elemental Microanalysis can be used, instead of the Leco® equivalent, each time a new combustion tube is required.

#### Benefits to the customer of the E4356 conversion kit

- Saves time because the lower O-ring needn't be removed during replacement.
- Allows more precise vertical alignment.
- Tube remains in lower tube ring during multiple cleanings, or until replaced.
- Easily installed, and can quickly be returned to original configuration.
- Considerable cost savings over the life of the CS844/744.



Our website features a video on how to use conversion kit E4356, – click the link to view the video or copy and paste the link into your browser:  
<http://elementalmicroanalysis.com/instructional-videos.php>



EMAL Part no	Description	OEM® part no
<a href="#">C4356</a>	Combustion tube CS744 and CS844 (Requires conversion kit E4356)	619-591-683
<a href="#">E4356</a>	Conversion kit for C4356 combustion tube (Contains 1x E4357, E4358, E4359, E1303, E1156)	
Each of the parts included in the conversion kit are also available separately:		
<a href="#">E4357</a>	Adaptor sleeve for C4356 combustion tube	
<a href="#">E4358</a>	Combustion tube ring insert for C4356 combustion tube	
<a href="#">E4359</a>	Insertion tool for C4356 combustion tube	
<a href="#">E1303</a>	O Ring	601-920
<a href="#">E1156</a>	O-Ring	

These products offer the same high quality and value for money that is typical of Elemental Microanalysis' entire product range. Available to order from our ISO 9001 factory online, by email, telephone or fax.

*All trademarks acknowledged. Use of the Leco® part number is for reference only and does not imply that the products are of Leco® origin. All Elemental Products are guaranteed to be of high quality and suitable for use in the stated application.*