

New product release

Leco 900/932 CHN, CHNS

Pre-Packed Combustion & Reduction Tubes

Save time and money by using our pre-packed tubes:

- ü Less downtime for your instrument
- ü No more measuring or handling of chemicals
- ü No waste
- ü Fast delivery
- ü High quality reagents



Tubes are supplied vacuum-packed to preserve reagent life-time.

Ordering details:

Cat No	Description	Pack
X3001	Prepacked Reduction tube Leco 900/932	1
X3002	Prepacked Reaction tube Leco 900/932 (CHN with Halogens)	1
X3003	Prepacked Reaction tube Leco 900/932 (CHN)	1
X3004	Prepacked Reaction tube Leco 932 (CHNS)	1
X3005	Prepacked Catalyst Heater tube Leco 900/932	1

Elemental Microanalysis
 1 Hameldown Road
 Okehampton
 Devon, EX20 1UB, UK
 +44 (0)1837 54446
info@microanalysis.co.uk
www.microanalysis.co.uk

Or your local dealer: