

Extended product range

EuroVector O rings

Elemental Microanalysis has extended its range of O rings and seals compatible with EuroVector analysers [EA3000](#), and [EA3100](#).

The additions to the product range enhance the existing product offerings for EuroVector including reagents, certified reference materials, glassware, sample containment, separation columns, prepacked tubes and other analysis accessories.

These O rings are identical to EuroVector products, Elemental Microanalysis Ltd can supply you from stock.

 E1511	 E1516	 E1519	 E1520
 E1522	 E1512	 E1510	 E1514
 E1523	 E1524	 E1517	 E1513
 E1515	 E1525	 E1518	 E1521

These products offer the same high quality and value for money that is typical of Elemental Microanalysis' entire product range.

Available to order from our [ISO 9001](#) factory online, by email, telephone or fax. Elemental Microanalysis offers extensive knowledge of analysis, over 40 years of experience and a worldwide reputation for excellence.

Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK

t. +44(0)1837 54446 • f. +44(0)1837 54544

www.elementalmicroanalysis.com • enquiries@microanalysis.co.uk



Available from your local Dealer:



Cross reference and analyser compatibility

EMAL part no.	New EuroVector part no.	Old EuroVector part no.	Product description	EA3000	EA3100
E1510	W08-031	EA213549	OR2050-silicon	X	X
E1511	W08-040	EA213554 EA213540-G	O ring 6mm for Q-fit	X	X
E1512	W08-041	EA213551	O ring for holder vectorSAS tightness(12,37 x 2,62mm) [since March 2013]	X	X
E1513	W08-042	EA213536	O ring 18mm neoprene for holder, Q-Fit & injector (since April 2014)	X	X
E1514	W08-043	EA213552	O ring for ash removal 13,10 x 2,62	X	
E1515	W08-044	EA213537	O ring 18mm for injector & ash removal	X	
E1516	W08-045	EA213542	EVR O ring for trap	X	
E1517	W08-111	EA213539	O ring 17mm for holder & Q-fit	X	
E1518	A02-009	EA213546	Parker seal 20mm	X	X
E1519	A04-054		Parker seal 8 mm		X
E1520	W08-029		8.75 x 1.78- ISO3601G		X
E1521	W08-035		26.5 x 2.65-ISO3601G		X
E1522	W08-036	EA213545	O ring for autosampler sample injector and viewing mirror tightness	X	X
E1523	W08-052		15.5 x 2.65-ISO3601G		X
E1524	W08-055	EA213534	17 x 1.8-ISO3601C		X
E1525	W08-102	EA213550	O ring for holder/vectorSAS tightness	X	

These products offer the same high quality and value for money that is typical of Elemental Microanalysis' entire product range.

Available to order from our [ISO 9001](#) factory online, by email, telephone or fax. Elemental Microanalysis offers extensive knowledge of analysis, over 40 years of experience and a worldwide reputation for excellence.

Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK
t. +44(0)1837 54446 • f. +44(0)1837 54544

www.elementalmicroanalysis.com • enquiries@microanalysis.co.uk



Available from your local Dealer:

