



## Product Sheet ~ Leco TruSpec Micro

### ~ Silver capsules/pans for solids

Cat #	Description	OEM #
D2300	Silver Capsules Pressed 4 x 3.2mm Light Weight pack of 500	502-206
D2000	Silver Capsules Pressed 4 x 3.2mm pack of 500	502-206

### ~ Tin capsules/pans for solids

Cat #	Description	OEM #
D1000	Tin Capsules Pressed Standard Weight 4 x 3.2mm pack of 500	502-227

### ~ Tools, trays etc

Cat #	Description	OEM #
E2017	Forceps Stainless Steel Straight Points	502-028
E2462	Spatula 501-614	501-614

### Chemical Reagents

Cat #	Description	OEM #
B1002	Copper Oxide Fine Wires 4 x 0.5mm 50gm	502-190
B1125	Copper Oxide Fine Wires 4 x 0.5mm 100gm	502-190
B1107	Copper Turnings Fine Long Strands Reduced and Annealed 501-621/502-656 200gm	501-621
		502-258
		502-295
		502-656
B1015M	Copper Wires Fine Wires Reduced 4 x 0.5mm 100gm	502-189
		502-304
		502-705
		502-844
B1212M	Copper Wires Fine Wires Reduced 4 x 0.5mm 1kg	-
B1083M	Copper Wires Pound Pack Fine Wires Reduced 4 x 0.5mm 454gm	502-304-500
B1161	EMADrone Magnesium Perchlorate 501-171 454gm	501-171
		501-171-HAZ
B1162	EMASorb A Granular 0.8 to 1.6mm 14 to 25 mesh 502-174 500gm	502-174
		502-174-HAZ
B1164	EMASorb B Granular 1.6 to 3mm 7 to 14 mesh 502-176 500gm	502-176
B1166	Glass Wool Borosilicate Fibres 501-081 454gm	501-081
B1167	Glass Wool Borosilicate Fibres 763-265 227gm	763-265
B1266	Glass Wool Rope Borosilicate Fibres 501-081 454g	501-081
B1267	Glass Wool Rope Borosilicate Fibres 763-265 227g	763-265
B1061	Nitrogen Catalyst Pellets 502-049 50gm	502-049
B1123	Rare Earth Copper Oxide Granular 0.85 to 1.7mm 501-170 50gm	501-170
B1188	Tungsten Trioxide Powder 502-301 125gm	502-301

### O-rings and operational accessories

Cat #	Description	OEM #
E2105	Absorbent Discs 3.5mm 502-283 pack of 100	502-283
E2154	Cleaning Brush 601-981	601-981
E2199	Lance Assembly Complete 616-147	616-147

E2197	Lance Holder Block 783-275	783-275
E2198	Lance Holder Cover 783-276	783-276
E1344	O Ring 601-684, 31.4mm x 2.6mm	601-684
E1345	O Ring 617-136, 5.2mm x 2.6mm	617-136
E1311	O Ring 771-407, 17.1mm x 2.6mm	771-407
E1342	O Ring 774-425, 41.0mm x 2.6mm	774-425
E1360	O ring 702-557, 29.8mm x 2.6mm	702-557
E1110	O Ring 772-520, 6.0mm x 2.6mm	772-520
E1101	O Ring Catalyst Heater Tube 601-442 14.0mm x 2.6mm	601-442
E1109	O Ring Doser Valve Assembly 772-269 pack of 10, 2.9mm x 1.8mm	772-269
E1114	O Ring Sample Drop Block 773-913, 4.5mm x 1.8mm	771-040 773-913
E1140	O-Ring Combustion Tube 616-138, 36.1mm x 3.5mm	616-138
E1139	O-Ring Lance Assembly 780-831 pack of 2, 17.2mm x 1.8mm	780-831
E2235	Secondary (screen) Filter 775-306	775-306
E2158	Silicone Grease 501-241 100gm	501-241
E2039	Snoop Leak Detection Fluid 236ml/8 fluid oz	502-213

### Pure Chemical Standards

Cat #	Description	OEM #
B2118	3,5-Dinitrobenzoic Acid OAS See Certificate no 308657 Expiry 27-Jun-23 5gm	502-204
B2114	Acetanilide OAS 5g See Certificate 293514 Expiry 15-Nov-22	501-053
B2104	Caffeine OAS 5gm See Certificate 239704 Expiry 30-april-24	502-205
B2115	Ethylenediaminetetra-acetic Acid (EDTA) OAS 502-092 50gm See Certificate 335390 Expiry 14-Jul-24	502-092
B2102	Stearic Acid OAS 5gm See Certificate 195037 Expiry 18-Dec-21	502-203
B2121	Sucrose OAS See Certificate 201037 Expiry 04-May-22 501-046 25gm	501-046
		501-441
B2106	Sulphamethazine OAS See Certificate 267171 Expiry 24-AUG-21 5gm	502-209

### Quartz, Glassware and Furnace tubes

Cat #	Description	OEM #
C4023	Catalyst Heater Tube Transparent Silica 601-437	601-437
C4338	Combustion Tube 619-065	619-065
C4358	Directing middle tube Truspec CHNS 619-777	619-777
C4027	Filling Funnel Borosilicate Glass 502-023	502-023
C4342	Filter tube 618-809	618-809
C4345	Glass tube TruSpec 619-425	619-425
C7000	Inner Pyrolysis Tube 608-130	608-130
C4366	Lance TruSpec CHNS Micro 619-774	619-774
C4051	Oxygen Lance Long (Lance length 5.25") 616-146	616-146
C4065	Porous Crucible (Pack of 10) 614-961-110	614-961 614-961-110
C4067	Porous Crucible (Pack of 5) (1/2 pack of 614-961-110)	-

## Product Sheet ~ Leco TruSpec Micro

		617-548	
C4332	Reagent Tube 617-548/633-103-225		633-103-225
C4337	Reagent Tube 6 774-638		774-638
C4341	Reagent tube 619-268		619-268
C4364	Reagent Tube 619-776		619-776
C4004	Reagent Tube Borosilicate Glass 775-601		775-601
C4339	Reduction Tube 619-154		619-154

**Quality:** All products are manufactured in our ISO9001:2015 approved facility and are covered by our 'Total Satisfaction Guarantee'.

**Can't see what you need?** We offer thousands of products for all elemental analysers. If you can't see what you need here the chances are we will have the product available. Just contact our friendly sales team who will track down what you need.

**OEM Part numbers:** The use of OEM part numbers is intended for convenience only and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and suitable for use in the stated instrument. All trade marks acknowledged.