


Products for Elementar® soli TOC cube



O-rings

Part No	Description	Pack	Elementar® part
E1054	Quadring, 7.65 mm x 1.78 mm	1	200005999, 05 000 353
E1054-10	Quadring, 7.65 mm x 1.78 mm	10	200005999, 05 000 353
E1067	O-ring, 4.5 mm x 1.5 mm	1	200008314, 05 000 249
E1067-10	O-ring, 4.5 mm x 1.5 mm	10	200008314, 05 000 249
E1069	O-ring, 7.5 mm x 2.5 mm	1	200006004, 05 000 416
E1069-10	O-ring, 7.5 mm x 2.5 mm	10	200006004, 05 000 416
E1087	O-ring, black, 7.59 mm x 2.62 mm	1	200014191, 05 001 302
E1171	O-ring for gripper arm	1	200010589, 03 002 395
E1172	O-ring, 32 mm x 3 mm	1	200008068, 03 002 268
E1173	O-ring, 5 mm x 1.5 mm	1	200006094, 05 000 425
E1175	O-ring, 16 mm x 3 mm	1	200007904, 05 000 568
E1175-10	O-ring, 16 mm x 3 mm	10	200007904, 05 000 568
E1555	Shaft seal	1	200006143, 03 002 607



Quartz & glassware

Part No	Description	Pack	Elementar® part
C3067	Combustion tube, quartz	1	200015419, 29.00-1078


Reagents & fillings

Part No	Description	Pack	Elementar® part
B1090	Alumina (aluminium oxide) spheres 3 to 6 mm	50 g	200010245, 50 008 467
B1354	Aquatak (drying agent) without indicator	50 g	200008689, 05 002 799
B1304	Brass wool	50 g	200010160, 38.00-0124
E2516	Cotton wool	50 g	200013169, 05 001 715
◆ B1141	Magnesium perchlorate, self indicating, granular, 0.7 to 1.2 mm	100 g	200006776, 05 001 059
B1255	Quartz wool discs, 24 mm diameter	25	
B1208	Silica (quartz) wool, very fine	10 g	200004827, 03 679 908


◆ These products are classified as hazardous for shipping



Brass wool
B1304



Magnesium perchlorate,
self indicating
B1141




The use of Elementar® part numbers is intended for convenience only, and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Products for Elementar® enviro TOC, vario TOC cube, vario TOC select


O-rings

Part No	Description	Pack	Elementar® part
E1054	Quadring, 7.65 mm x 1.78 mm	1	200005999, 05 000 353
E1054-10	Quadring, 7.65 mm x 1.78 mm	10	200005999, 05 000 353
E1067	O-ring, 4.5 mm x 1.5 mm	1	200008314, 05 000 249
E1067-10	O-ring, 4.5 mm x 1.5 mm	10	200008314, 05 000 249
E1069	O-ring, 7.5 mm x 2.5 mm	1	200006004, 05 000 416
E1069-10	O-ring, 7.5 mm x 2.5 mm	10	200006004, 05 000 416
E1072	O-ring, 18 mm x 2 mm	1	200013028, 05 000 096
E1073	O-ring, rubber, 20 mm x 2 mm	1	200007224, 03 654 627
E1073-10	O-ring, rubber, 20 mm x 2 mm	10	200007224, 03 654 627
E1075	O-ring, rubber, 24 mm x 2 mm	1	200007469, 03 002 267
E1075-10	O-ring, rubber, 24 mm x 2 mm	10	200007469, 03 002 267
E1087	O-ring, black, 7.59 mm x 2.62 mm	1	200014191, 05 001 302
E1172	O-ring, 32 mm x 3 mm	1	200008068, 03 002 268
E1556	O-ring, 18 mm x 2.5 mm	1	200004884, 05 001 053



Quartz & glassware

Part No	Description	Pack	Elementar® part
C3054	Ash crucible, quartz, 60 mm	1	200014796, 38.00-1290
C3018	Ash finger, quartz, 60 mm	1	200014555, 11.00-1170/4
C3060	Ash finger, quartz, 60 mm, pack of 10	10	200014555, 11.00-1170/4
C3009	Combustion/reduction tube	1	200007844, 11.00-1150/4
C1143	Drying tube 300 mm (tube only, no caps)	1	200008714, 15.00-0080
C3026	Drying tube 348 mm (with GL25 caps)	1	200011489, 16.00-1012
C3045	Filter tube	1	200013470, 50 009 375/4
C3074	Filter tube 42 mm	1	200007012, 16.00-1081
C3021	Protection tube 105 mm	1	200014026, 11.00-1317/4
C3059	Protection tube 105 mm (alt to C3021 without slits)	1	200014026, 11.00-1317/4



The use of Elementar® part numbers is intended for convenience only, and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Products for Elementar® enviro TOC, vario TOC cube, vario TOC select (continued)


Reagents

Part No	Description	Pack	Elementar® part
B1089	Alumina (aluminium oxide) wool, fine	20 g	200004773, 50 000 086
B1354	Aquatak (drying agent) without indicator	50 g	200008689, 05 002 799
B1304	Brass wool	50 g	200010160, 38.00-0124
B1270	Cerium dioxide, fused, granular	100 g	200013646, 05 000 469
B1218	Copper oxide wires (coarse), 6 x 0.65 mm	250 g	200007554, 03 679 905
B1122	Copper oxide wires (fine), 4 x 0.5 mm	250 g	200009554, 05 001 039
◆ B1141	Magnesium perchlorate, self indicating, granular, 0.7 to 1.2 mm	100 g	200006776, 05 001 059
B1196	Platinum catalyst, 5% Pt	30 g	200007811, 03 002 262
B1308	Platinum catalyst, 5% Pt	150 g	200007811, 03 002 262
B1601	Potassium hydrogen phthalate	25 g	200009208, 35.00-0151
B1255	Quartz wool discs, 24 mm diameter	25	
B1215	Silica (quartz) chips, granular, 2 to 4 mm	50 g	200007521, 05 000 900
B1208	Silica (quartz) wool, very fine	10 g	200004827, 03 679 908
B1231	Silver wool, 0.05 mm fine wire	50 g	200005872, 22 131 365
B1602	Sodium carbonate	25 g	200007465, 35.00-0152

Brass wool
B1304



Cerium dioxide, fused, granular
B1270



Platinum catalyst, 5% Pt
B1196, B1308



◆ These products are classified as hazardous for shipping

Operating supplies

Part No	Description	Pack	Elementar® part
E1174	Filter, diameter 18 mm	1	200015256, 03 002 399
E1174-10	Filter, diameter 18 mm	10	200015256, 03 002 399
E2533	Set of ball valve half shells 14 mm	1	200010811, 21.00-0008

Filter, diameter 18 mm
E1174



Set of ball valve half shells, 14 mm
E2533



Accessories

Part No	Description	Pack	Elementar® part
E2516	Cotton wool	50 g	200013169, 05 001 715
E2016	Forceps, stainless steel, curved points - length 130 mm	1	200010676, 03 001 248
E2538	GL14 blanking cap, red	2	
E2518	GL18 screw cap	1	200011744, 05 000 044
E2505	Gloves	Pair	200006785, 05 000 094
E2002	Macro multi-purpose spoon/spatula 215 mm	1	
E2521	Manually operated drill	1	200009657, 11.00-0500/4
E2064	Multiwell plate, 24 cell, flat bottom	1	200013850, 05 000 308
E2526	Powder filling funnel (small) with silicone tubing outlet	1	200004594, 34.00-0051/4
E2240	Vacuum grease	50 g	200006998, 03 968 006

GL18 screw cap
E2518



GL14 blanking cap
E2538





The use of Elementar® part numbers is intended for convenience only, and does not imply that the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Products for Elementar® vario TOC cube, vario TOC select (continued)

Sample encapsulation

Part No	Description	Pack	Elementar® part
D2008	Silver capsules pressed 8 x 5 mm	250	200014521, 05 001 540
D4103	Silver capsules smooth wall 5 x 2 mm	250	200013710, 24 005 100
D4102	Silver capsules smooth wall 5.5 x 3.5 mm	250	200007516, 05 000 816
D4109	Silver capsules smooth wall 9 x 3.5 mm	250	200010387, 05 000 817
D5127	Silver weighing boats 6 x 4 x 4 mm	500	200008833, 11.02-1046
D5039	Silver weighing boats 12 x 4 x 4 mm	500	200013336, 11.00-1371
D5026	Silver weighing boats 12 x 6 x 6 mm	100	200007754, 22 133 212
D5034	Silver weighing boats 15 x 8 x 8 mm	70	200013395, 22 133 213
D5035	Silver weighing boats 20 x 10 x 9 mm	42	200010913, 22 133 214
D4033	Smooth wall tin capsules flat base 5.5 x 3.5 mm	100	200013224, 05 001 727
D4131	Smooth wall tin capsules flat base 6 x 3 mm	250	200014058, 03 951 620
D4132	Smooth wall tin capsules flat base 9 x 3.5 mm	250	200015479, 03 951 619
D4121	Smooth wall tin capsules flat base 11 x 4.75 mm	250	200008184, 05 000 481
D4120	Smooth wall tin capsules flat base 14 x 6.75 mm	100	200013285, 05 000 312
D4009	Smooth wall tin capsules round base 20 x 8 mm	100	200011997, 05 000 311
D5029	Tin boats 6 x 4 x 4 mm	1000	200007543, 11.02-1047
D5130	Tin boats 12 x 4 x 4 mm	1000	200014577, 22 137 418
D5131	Tin boats 12 x 6 x 6 mm	500	200007074, 22 137 419
D5032	Tin boats 15 x 8 x 8 mm	250	200012350, 22 137 420
D5033	Tin boats 20 x 10 x 9 mm	120	200006296, 12.00-1009
D5140	Tin boats heavy weight 12 x 4 x 4 mm	1000	200007257, 11.00-5057
D9001	Tin foil cones standard weight	100	200009010, 05 000 429
D1068	Tin foil squares standard weight 37 x 37 mm	100	200009010, 05 000 429
D1118	Tin foil squares standard weight 50 x 50 mm	500	200011221, 12.01-0033
D1119	Tin foil squares standard weight 50 x 50 mm	1500	200008484, 34.00-0052



The use of Elementar® part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Elementar® trademarks acknowledged.

Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK
 • t. +44(0)1837 54446 • www.elementalmicroanalysis.com
 • enquiries@microanalysis.co.uk



Available from your local Dealer:

