

## Products for Antek MultiTek<sup>®</sup>, 9000 Series & 7000 Series


### Quartz & glassware

| Part No | Description                  | Pack | Antek part |
|---------|------------------------------|------|------------|
| C3110   | Pyrotube, direct inject      | 1    | 70061      |
| C3111   | Pyrotube, dual pass 9000 LLS | 1    | 71303      |
| C3112   | Pyrotube, multi-matrix       | 1    | 71121      |
| C3113   | Pyrotube, multi-matrix       | 1    | 71127      |
| C3122   | Pyrotube, 90° elbow          | 1    | 70008      |
| C3115   | Quartz chips (pack = 1 fill) | 67 g | 71029      |



### Quartz boats for 9000/7000

| Part No | Description                             | Pack | Antek part |
|---------|---|------|------------|
| C3117   | Sample boat, liquids, multi-matrix      | 1    | 71122      |
| C3118   | Sample boat, water jacket, Pyro         | 1    | 70007      |
| C3116   | Sample weigh boat carrier, multi-matrix | 1    | 71124      |
| C3119   | Sample weigh boat, multi matrix         | 1    | 71123      |
| C3121   | Quartz weigh boat, solids               | 1    | 71039      |




## Antek MultiTek<sup>®</sup>, 9000 Series & 7000 Series

### Other consumables & parts

| Part No | Description                         | Pack  | Antek part |
|---------|-------------------------------------|-------|------------|
| E7902   | Cera-wool 2" x 2"                   | 1     | 79127      |
| B1715   | Drierite, H <sub>2</sub> O removing | 1 lb  | 71028      |
| E7901   | Filter, parker                      | 1     | 71111      |
| E7910   | Membrane dryer                      | 1     | 70058      |
| B1710   | Moleculite, O <sub>3</sub> scrubber | 1 b   | 72001      |
| E1560   | O-ring, tube inlet                  | 1     | 71194      |
| E1561   | O-ring, tube exit                   | 1     | 70256      |
| E7900   | Septum, 12.5 mm                     | 10    | 70002      |
| E7905   | Septum, 12.5 mm                     | 100   | 69091      |
| E7906   | Syringe, 100 µl                     | 1     | 90073      |
| E7907   | Syringe, 25 µl, conetip             | 1     | 36117      |
| E7908   | Teflon ferrules, ¼"                 | 3     | 70027      |
| E7909   | Teflon tubing, ¼"                   | 50 ft | 30066      |



Moleculite  
B1710



Drierite  
B1715

Cera-wool  
E7902



Teflon ferrules  
E7908



E1560



E1561

Septum, 12.5 mm  
E7900, E7905



Syringe, 100 µl  
E7906



Parker filter  
E7901



Syringe, 25 µl, conetip  
E7907



Membrane dryer  
E7910



The use of Antek part numbers is intended for convenience only, and does not imply the products are of OEM origin. All Elemental Microanalysis products are guaranteed to be of high quality and are deemed to be suitable alternatives in the stated instruments. All Antek trademarks acknowledged.

### Elemental Microanalysis Ltd

1 Hameldown Road, Okehampton, EX20 1UB, UK  
 • t. +44(0)1837 54446 • [www.elementalmicroanalysis.com](http://www.elementalmicroanalysis.com)  
 • [enquiries@microanalysis.co.uk](mailto:enquiries@microanalysis.co.uk)



Available from your local Dealer:

